

Product Range

General Purpose Dry Wipes

Product Type

Hygenex Eco-Friendly Patient Dry Wipes

Distributor

DDC Dolphin Ltd
 The Fulcrum
 Vantage Way
 Poole
 Dorset
 BH12 4NU
 +44 (0)1202 731555
 info@hygenex.co.uk
 www.hygenex.co.uk



Product Part Code

90-01-088

Country of Manufacture

United Kingdom

Unit of Sale

Product Number	90-01-088
Unit of Sale	Pack / Case / Pallet
Pack Quantity	75 Wipes
Pack Dimensions (mm)	260 x 140 x 60
Pack Weight (g)	340
EAN Pack Barcode	5060542000413
Case Quantity	16 Packs
Case Dimension (mm)	580 x 320 x 320
Case Weight (kg)	6
EAN Case Barcode	5060542000451
Wipe Area (cm ² ± 10%)	674
Pallet Dimensions (mm)	1200 x 1000 x 1635
Total Loaded Pallet Weight (kg)	200
Pallet Height (mm)	1635
EAN Pallet Barcode	5060542000475
Customs Stat No (HS Code)	56031210

HYGENEX ECO-FRIENDLY PATIENT DRY WIPES

Product Information

Product Description

Dry patient and general-purpose wipes that have been specially developed to deliver outstanding performance whilst offering excellent environmentally friendly benefits. Soft when dry and yet strong and smooth when wet, making them kind to patients' skin and highly effective in their ability to clean. Highly absorbent, the wipes still retain strength and integrity in use. 100% biodegradable, flushable, dispersible and maceratable

Product Features

- 100% biodegradable, compostable, dispersible, flushable and maceratable.
- Safe to dispose in toilets, macerators, bedpan washers or slop hoppers.
- Dual application wipes for use on surfaces or skin.
- Suitable for use with detergents, disinfectants, and skin cleansing products.
- Kind to skin and highly effective in their ability to clean.
- Highly absorbent.
- Suitable for full body cleansing (Individual wipe size 270 x 260mm).

Intended Uses

- Intended for general surface and skin cleansing.
- Ideal for lifting urine, faeces and menses from soiled skin and surfaces.
- Suitable for use with detergents, disinfectants, and skin cleansing products.
- Suitable for professional use in hospitals, nursing homes, care home and domiciliary care and also for personal care in the home or at work.
- Suitable for use on or by patients and also for occupational use by care workers to cleanse their skin, wipes are soft when dry and smooth when wet.

Non-woven Material

- 70 gsm 100% cellulosic material.
- Disposal - 100% biodegradable, flushable, dispersible, maceratable.
- Low to mid strength material. Generally, non-linting but can lint on facial stubble.
- Product disperses following normal agitation in water during flushing and drain transit. It can be macerated in all macerators and will clear domestic and municipal pumps. It biodegrades under aerobic and anaerobic conditions.
- Cellulosic content derived from sustainable forestry. FSC and PEFC approved.
- Care is required on dispensing due to construction from 100% percent naturally derived material.
- Product meets EDANA/INDA standards for flush ability which require that the following tests are passed.
- Toilet and drain line clearance test (will it flush and transit along drain line).
- Slop box disintegration test (dispersion).
- Household pump test (anti clog).
- Settling test (anti float).
- Aerobic bio disintegration and bio degradation test.
- Anaerobic bio disintegration and bio degradation test.
- Municipal sewage pump test (anti clog).

Wipes Presentation, Dispensing & Use

- Peel back the perforated opening to access the wipes.
- Wipes are 'Z' folded and interleaved for single dispensing and additionally half folded for packing efficiency.
- Wipe seam is presented in the middle of the aperture for ease of dispensing.
- Wipe seam is always the same direction to aid dispensing. With the pack the correct way around to read the information on it, stroke towards the user with fingertip to elevate edge of wipe. Then pinch lightly between fingertips and gently remove from pack.

HYGENEX ECO-FRIENDLY PATIENT DRY WIPES

- Product Information

- Take care with dispensing cellulosic wipes (flushable/dispersible/maceratable types). The natural materials used in their construction cause a slight adherence between wipes and aggressive removal will cause ripping.
- Take care with removal of the first few wipes as they are packed most tightly in the pack.
- Dispose of responsibly.
- Use each wipe once only.



DDC Dolphin Ltd

The Fulcrum, Vantage Way, Poole, Dorset BH12 4NU

tel: **01202 731555** | fax: 01202 731666 | info@ddcdolphin.com

DDC Dolphin and Hygenex are Registered Trade Marks