

# MIDI SLIPPER PAN LINER

# HYGENEX

## DISPOSABLE MIDI SLIPPER (FRACTURE) PAN LINER 1.3LTR



### PRODUCT DESCRIPTION

Disposable Slipper Pan Liner, made of pulp material with water resistance properties.

1.3 litre capacity.

Suitable for disposal in a medical plug macerator.

### INTENDED USE

The disposable slipper pan is used for holding human waste and other bodily fluids. It is also used by immobile female patients for bed-side collection of urine.

Designed to be used with a slipper pan support to easily slide under a patient and accept body weight up to 160kgs (DDC Dolphin Part Number 91-02-009).

### KEY BENEFITS

- An environmentally friendly product that is made of recycled newspaper
- Single use and disposable – significantly reducing the risk of cross contamination and Healthcare Associated Infections
- Water resistant for a minimum of 4 hours
- Stackable for storage purposes and easy separation for dispensing
- Environmentally friendly cardboard packaging
- Can be weighed to record fluid output

Product code 90-01-003

Contact us today:



+44(0) 1202 731555



[sales@hygenex.co.uk](mailto:sales@hygenex.co.uk)



[www.hygenex.co.uk](http://www.hygenex.co.uk)

# MIDI SLIPPER PAN LINER

# HYGENEX

## PRODUCT SPECIFICATION

Manufacturer's Product Code	91-01-003
Unit of Issue – Carton	100 units per carton
Unit of Issue – Pallet	17 cartons per pallet
Usable Capacity	1300ml
Product Dimensions - L x W x H (Approx.)	350 x 285 x 75mm
Product Weight	46 ± 4g
Carton Dimensions – L x W x H (Approx.)	355 x 300 x 805mm
Carton Weight (Approx.)	5.1kg
Pallet Dimensions – L x W x H (Approx.)	1200 x 1050 x 1250mm
Pallet Weight (Approx.)	92kg
EAN 13 Barcode	5060542000031

## PRODUCT INFORMATION

Composition	Clean, over issued newspaper with fixative, waterproofing additives
Shelf Life	Minimum 3 years if stored in a dry environment with a temperature range of -5°C to +40°C
Sterilisation	Non-sterile
Classification	Class 1 (According to Annex 1X, Rule 1 of Council Directive 93/42EEC, amended by directive 2007/47/EC)
Standards and Guidelines	PAS29:1999 Disposable pulp products for use in healthcare (maceration test normal cycle) PAS29:199 Disposable pulp products for use in healthcare, clause 5 (water retention test at 35 degrees centigrade +/- 3%) ISO9001:2015 BRC / IOP Global Standard
Disposal Method	Macerator or in accordance with local waste policy
Country of Origin	United Kingdom
Harmonized Tariff Commodity Code	4823 7090

Contact us today:



+44(0) 1202 731555



[sales@hygenex.co.uk](mailto:sales@hygenex.co.uk)



[www.hygenex.co.uk](http://www.hygenex.co.uk)